

Manual Crimping Tools



CRS's ergonomically curved handle design conforms to your natural hand grip, which greatly improves comfort over blocky metal grips. The body is constructed of tough, light-weight glass-filled plastic for strength.

The bottom-pull handle design is steadier to hold than top-push handle designs. The head does not wobble. It is easy to hold the tapered jaws of the crimper in a stable position, as the bottom handle is pulled up.

The adjustment knob is conveniently positioned at the head of the crimper, clearly displays + and - symbols along with directional arrows. This simplifies adjusting the intensity of the desired crimp. When the knob touches the crimper body, the crimp setting has been reached. No more "extra squeeze" necessary!

To ensure a smooth and consistent operation, the hand crimper uses the same precision tolerances as the CRS Electronic Crimper and Decapper.

Description	Item No.
11 mm Manual Crimper	311990
20 mm Manual Crimper	320990
13 mm Manual Crimper	313990
8 mm Manual Crimper	308990

Description	Item No.
11 mm Manual Decapper	311991
20 mm Manual Decapper	320991
13 mm Manual Decapper	313991
Base	5ASTA0

20 mm Manual Crimper, for Flip Off Caps	320992
13 mm Manual Crimper, for Flip Off Caps	313992